












Side cutting spatula double armed			Metric	USP	Length	Order number
1/2 Circle 	2 x HSP 8 mm 200 µm		0.5	7/0	45 cm	<i>violet</i> OSU 07
Side cutting spatula double armed			Metric	USP	Length	Order number
3/8 Circle 	2 x DSP 6 mm 200 µm		0.4	8/0	45 cm	<i>violet</i> OSU 08
	2 x DSP 8 mm 200 µm		0.5	7/0	45 cm	<i>violet</i> OSU 071
Side cutting spatula double armed			Metric	USP	Length	Order number
1/4 Circle 	2 x VSP 8 mm 300 µm		0.7	6/0	45 cm	<i>violet</i> OSU 06
	350 µm		1	5/0	45 cm	<i>violet</i> OSU 05