



MICROSURGERY











*For quality
&
safety*










www.sutures.be

Microsurgery
Revision 4

MICRO-
SURGERY

Micro needle taper point 	Metric	USP	Length	Order number			
3/8 Circle 	MDR 3,3 mm 50µm		0.1	11/0	10 cm	<i>black</i>	MN 111
	MDR 3,8 mm 50µm		0.2	10/0	10 cm	<i>black</i>	MN 102
	MDR 5,1 mm 70µm		0.01	12/0	10 cm	<i>black</i>	MN 124
			0.1	11/0	10 cm	<i>black</i>	MN 114
			0.2	10/0	10 cm	<i>black</i>	MN 104
	MDR 3,8 mm 70µm		0.1	11/0	10 cm	<i>black</i>	MN 113
	MDR 3,8 mm 100µm		0.2	10/0	15 cm	<i>black</i>	MN 105
	MDR 5,1 mm 100µm		0.2	10/0	15 cm	<i>black</i>	MN 106
			0.3	9/0	15 cm	<i>black</i>	MN 096
	MDR 6,4 mm 140µm		0.4	8/0	15 cm	<i>black</i>	MN 087
MDR 9,7 mm 200µm		0.5	7/0	45 cm	<i>black</i>	MN 078	

Micro needle taper point		Metric	USP	Length	Order number	
3/8 Circle 	 MDR 5,1 mm 100µm		0.2	10/0	15 cm	<i>blue</i> MP 106
	MDR 6,4 mm 150µm		0.3	9/0	15 cm	<i>blue</i> MP 097

Micro needle taper point		Metric	USP	Length	Order number
3/8 Circle 	 MDR 6,4 mm 150µm	 0.4	8/0	20 cm	<i>violet</i> M910 087